

# Raystation Scripting

Brian Mark Anderson

4/2/2020

RayStation 7

Patient name: BMA\_Example ID: 000000011 Case: CASE 1

RayStation 6.99 Not For Clinical Use - Evaluation Use Only

Automated Planning

Patient Data

Patient Modeling

Plan Design

Optimization

Evaluation

QA

Adaptation

Patient Information

Open

Import new patient

Delete

Save

Close

Edit patient data

Edit case data

Import

Export

Create plan report

Initialize treatment adaptation

PATIENT HANDLING

CURRENT

ROI

Organs at risk (1)

External

PATIENT INFORMATION

CURRENT PLAN

TREATMENT ADAPTATION

CURRENT CASE

IMAGE SETS

Properties

Delete

CT 1

CT 2

CT 3


Create new material...

Remove material...



# First things first

Console can do everything





StateTree is your best friend.. Everything that can be done On Raystation is in the state tree

```
IronPython 2.7.1 (2.7.0.40) on .NET 4.0.30319.42000
Type "help", "copyright", "credits" or "license" for more information.
Connecting to RayStation. (Session id = 11080_2 )
Variable 'case' is set to case with name: "CASE 1"
Variable 'examination' is set to examination with name: "T1"
Variable 'machine_db' is set to the current machine database
Variable 'patient_db' is set to the current patient database
Variable 'ui' is set to the current ui
Imported statetree
>>> statetree.RunStateTree()
```

- Cases
- ImportLogs
- ModificationInfo
  - CalculateGammaForFractionDose(...)
  - EditPatientInformation(...)
  - EditRoiLineWidthVisualization(...)
  - EditShowDoseGridVisualization(...)
  - EditShowDoseVisualization(...)
  - GetImagingSystemsNames(...)
  - ImportDicomDataFromPath(...)
  - ImportDicomDataFromRepository(...)
  - Save(...)
  - Set2DVisualizationForRoi(...)
  - SetPoiVisibility(...)
  - SetRoiVisibility(...)
  - SetShowDoseGridOutlineVisualization(...)
  - SetShowDoseGridVoxelsVisualization(...)

# Lets make an ROI





ROI Tools

POI Tools

Approval

Fusion

ROI

Matl

▼ Organs at risk (1)

Test\_1

▼ Unknown (2)

Test

Test2

RayStation Command Console -X:TabCompletion -X:ColorfulConsole -X:AutoIndent -X:Frames




File Edit View Help

IronPython 2.7.1 (2.7.0.40) on .NET 4.0.30319.42000  
Type "help", "copyright", "credits" or "license" for more information.  
Connecting to RayStation. (Session id = 11892\_1 )  
Variable 'case' is set to case with name: "CASE 1"  
Variable 'examination' is set to examination with name: "CT 1"  
Variable 'machine\_db' is set to the current machine database  
Variable 'patient\_db' is set to the current patient database  
Variable 'ui' is set to the current ui  
Imported statetree  
>>> statetree.RunStateTree()  
>>> case.PatientModel.CreateRoi(Name='Test')  
<ScriptObject id=7f1395a4-b3d7-41d3-8b73-044de8ca97bb, 'Test'>  
>>> case.PatientModel.CreateRoi(Name='Test',Color='Blue')  
Traceback (most recent call last):  
File "<stdin>", line 1, in <module>  
SystemError: Name 'Test' is not unique  
>>> case.PatientModel.CreateRoi(Name='Test2',Color='Blue')  
<ScriptObject id=5c5c2110-870d-4bc1-9f19-b81692e65e94, 'Test2'>  
>>>

Ready 25x80

# How to know what to type?

- Look in statetree, or cross fingers and hope you can record it



# Results of Record

This can at least give you a hint of where to look in the statetree. Note that this is case.PatientModel

```
# Script recorded 21 Jun 2018

# RayStation version: 6.99.0.14
# Selected patient: ...

from connect import *

case = get_current("Case")

retval_0 = case.PatientModel.CreateRoi(Name="Temp3", Color="Red", Type="Organ", TissueName=None, RbeCel:
```

# Making A Script


```
1  from connect import *
2
3  case = get_current('Case')
4  case.PatientModel.CreateRoi(Name="Test3", Color="Blue", Type="Organ",
5                               TissueName=None, RbeCellTypeName=None,
6                               RoiMaterial=None)
```



# Making External ROI

- Can be helpful for rigid registrations
- Record action\* But we want to repeat it across all exams

# Make External ROI

```
1  from connect import *
2
3  case = get_current('Case')
4  patient = get_current('Patient')
5
6  case.PatientModel.CreateRoi(Name="External", Color="Green",
7                               Type="External", TissueName="",
8                               RbeCellTypeName=None, RoiMaterial=None)
9
10 case.PatientModel.RegionsOfInterest['External'].\
11      CreateExternalGeometry(Examination=case.Examinations['CT 1'], ThresholdLevel=-250)
```

# Make External ROI across exams

```
1  from connect import *
2
3  case = get_current('Case')
4  patient = get_current('Patient')
5
6  case.PatientModel.CreateRoi(Name="External", Color="Green",
7                               Type="External", TissueName="",
8                               RbeCellTypeName=None, RoiMaterial=None)
9
10 for exam in case.Examinations:
11     case.PatientModel.RegionsOfInterest['External'].\
12         CreateExternalGeometry(Examination=exam, ThresholdLevel=-250)
```

# Make External ROI across Cases

```
1  from connect import *
2
3  case = get_current('Case')
4  patient = get_current('Patient')
5
6  for case in patient.Cases:
7      case.PatientModel.CreateRoi(Name="External", Color="Green",
8                                  Type="External", TissueName="",
9                                  RbeCellTypeName=None, RoiMaterial=None)
10     for exam in case.Examinations:
11         case.PatientModel.RegionsOfInterest['External'].\
12             CreateExternalGeometry(Examination=exam, ThresholdLevel=-250)
```

# Make External across Patients and cases

- First need to load up a new patient..
- Define a way to query patients:

```
from connect import *  
  
def ChangePatient(patient_db, patient_id):  
    info = patient_db.QueryPatientInfo(Filter={"PatientID": patient_id}, UseIndexService=True)  
    patient = patient_db.LoadPatient(PatientInfo=info[0], AllowPatientUpgrade=True)  
    return patient
```

# Load a patient

```
1  from connect import *
2
3
4  def ChangePatient(patient_db, patient_id):
5      info = patient_db.QueryPatientInfo(Filter={"PatientID": patient_id}, UseIndexService=True)
6      patient = patient_db.LoadPatient(PatientInfo=info[0], AllowPatientUpgrade=True)
7      return patient
8
9  patient_db = get_current("PatientDB")
10 MRNs = ['000000011']
11 for MRN in MRNs:
12     patient = ChangePatient(patient_db, MRN)
13     for case in patient.Cases:
14         case.SetCurrent()
```

```

from connect import *

def ChangePatient(patient_db, patient_id):
    info = patient_db.QueryPatientInfo(Filter={"PatientID": patient_id}, UseIndexService=True)
    patient = patient_db.LoadPatient(PatientInfo=info[0], AllowPatientUpgrade=True)
    return patient


def main():
    patient_db = get_current("PatientDB")
    MRNs = ['000000011']
    case = get_current('Case')
    for MRN in MRNs:
        patient = ChangePatient(patient_db, MRN)
        for case in patient.Cases:
            case.PatientModel.CreateRoi(Name="External", Color="Blue", Type="External",
                                         TissueName=None, RbeCellTypeName=None,
                                         RoiMaterial=None)
            for exam in case.Examinations:
                case.PatientModel.RegionsOfInterest._External.CreateExternalGeometry(
                    Examination=exam)

if __name__ == '__main__':
    main()

```

# Script create a new ROI

- Open the console
- Run `statetree.RunStateTree()`
- See all available options
- Remember, we were under the `PatientModel` in raystation







# Scriptability of actions

- case.ComputeRigidImageRegistration(FloatingExaminationName=Reference\_Image,  
ReferenceExaminationName=exam.Name,  
UseOnlyTranslations=False,  
HighWeightOnBones=True,  
InitializeImages=True,  
FocusRoisNames=[],  
RegistrationName="**BMA\_Rigid**")

Thank you!